

# How to set up a lab balance in 6 steps

Proper balance setup is imperative to help ensure consistent, accurate results. Here are 6 steps to help you set up your balance correctly.

# 1

Always read the product manual before setting up the balance.

# 2

Place the balance on a solid, stable surface, away from direct sunlight and any sources of vibration, breezes and moisture.

# 3

Verify that the balance is sitting evenly on the surface by checking the balance's leveling bubble.

# 4

Adjust the balance feet by raising or lowering them one at a time until the bubble is centered.

# 5

An anti-vibration table can help increase stability of the balance, improving accuracy.

# 6

For the best results, make certain the room temperature remains consistent when the balance is in use.

