
MTB SCALE (EN)
BALANCE MTB (FR)
MTB WAAGE (DE)
MTB BILANCIA DIGITALE (ES)
SERIE MTB (IT)

(P.N. 7.00.6.6.0263 - Revision A1, September 2013)

ENGLISH:	P	1	-	3
FRANÇAIS:	P	4	-	6
DEUTSCH:	P	7	-	9
ESPAÑOL:	P	10	-	11
ITALIANO:	P	12	-	15
EU DECLARATION:	P	16		

1. Installation diagram

SPECIFICATIONS

Capacity (max.)	20 kg 44 lb
Readability	0.005 kg 0.005 lb
Weighing units	kg and lb
Power	6V, 500mA DC adapter
Batteries	4 x AA size
Pan Size	560 x 285 mm (22" x 11.2")
Overall Size (w x h x d)	560 x 350 x 175 mm (22" x 13.75" x 7.0")
Net weight	3.2 kg (7 lb.)

OPERATION

Install batteries, ensuring the polarity is as shown in the battery holder, or insert the power adapter jack into the connector on the rear of the scale and turn on.

Place the scale on a flat, very stable surface, and clear the scale platform.
Press the **[ON/ZERO]** key to turn on the scale.

Display will show 0.000 kg or 0.000 lb. depending on the last weighing unit selected.
It is now ready to weigh.

Put the subject on the scale platform, as central as possible. The display shows the weight.
Once reading is stable, the display flashes.
Press the **[ON/ZERO]** key briefly to tare the scale, setting the display to 0.000.
Any additional weight now added will displayed.

SELECT WEIGHING UNIT

Press the **[MODE]** key to change from kilograms to pounds and back again.
The display will show the weighing unit selected.

STORE THE WEIGHT

To store the displayed weight: Press the **[MEM]** key when the LCD display is flashing.

RECALL THE STORED THE WEIGHT

To recall the stored weight: Press the **[MEM]** key when the LCD display is at zero
Press it a second time to return to weighing.

To turn off the scale.

Press and hold the **[ON/ZERO]** key

The scale will turn off automatically after 5 minutes if it is not being used.

CALIBRATION

Place the scale on a flat, very stable surface, and clear the scale platform.

When the scale is off press the **[MODE]** key, then press the **[ON/ZERO]** key. Hold both keys down until the display shows "CAL" then release both keys.

The display will show a number representing the internal count from the sensor
Press the **[ON/ZERO]** key briefly.

The display will flash "**C10**" and then show " 0".

Place an accurate reference 10kg mass on the platform (OIML class M1 recommended)

The display will show a number from the sensor.

Press the **[ON/ZERO]** key.

The display will flash "**C20**" and then show " 0000".

Press the **[ON/ZERO]** key again, the display will show “End” then turn off.
The calibration is complete.

Remove the 10kg mass from the platform and turn on the scale.
Place the 10kg on the scale, it should display 10.000kg $\pm$ 0.010kg. If it is not within acceptable values repeat the calibration.

WARNING MESSAGES

“**Loo**”: Low power, please change the batteries.

“**o_Ld**”: Overloading (>23 kg, 50 lb), please remove the item immediately to protect the scale.

WARRANTY STATEMENT

Adam Equipment offers Limited Warranty (Parts and Labour) for any components that fail due to defects in materials or workmanship. Warranty starts from the date of delivery.

During the warranty period, should any repairs be necessary, the purchaser must inform its supplier or Adam Equipment Company. The company or its authorised Technician reserves the right to repair or replace the components at any of its workshops at no additional cost, depending on the severity of the problems. However, any freight involved in sending the faulty units or parts to the Service Centre should be borne by the purchaser.

The warranty will cease to operate if the equipment is not returned in the original packaging and with correct documentation for a claim to be processed. All claims are at the sole discretion of Adam Equipment.

This warranty does not cover equipment where defects or poor performance is due to misuse, accidental damage, exposure to radioactive or corrosive materials, negligence, faulty installation, unauthorised modifications or attempted repair, or failure to observe the requirements and recommendations as given in this User Manual.

This product may include a rechargeable battery that is designed to be removed and replaced by the user. Adam Equipment warrants that it will provide a replacement battery if the battery manifests a defect in materials or workmanship during the initial period of use of the product in which the battery is installed.

As with all batteries, the maximum capacity of any battery included in the product will decrease with time or use, and battery cycle life will vary depending on product model, configuration, features, use, and power management settings. A decrease in maximum battery capacity or battery cycle life is not a defect in materials or workmanship, and is not covered by this Limited Warranty.

Repairs carried out under the warranty do not extend the warranty period. Components removed during warranty repairs become company property.

The statutory rights of the purchaser are not affected by this warranty. The terms of this warranty is governed by the UK law. For complete details on Warranty Information, see the terms and conditions of sale available on our web-site.

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16

Year CE marking was first affixed to declared product

MTB DIGITAL ELECTRONIC WEIGHING SCALE SERIES

Maidstone Road, Kingston
Milton Keynes, MK10 0BD

United Kingdom

Adam Equipment Co. Ltd.

Tel: +44 (0) 1908 274545

Fax: +44 (0) 1908 641339

Email: sales@adamequipment.co.uk

www.adamequipment.com

We declare under our sole responsibility that the balance models listed below marked with "CE" are in conformity with the directives and standards mentioned.

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Signed for and on behalf of: Adam Equipment Company Ltd.

United Kingdom, 20 Feb 2016

Clive Jones, Quality & Product Manager.

WEEE 2012/19/EU

This device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements. Disposal of batteries (if fitted) must conform to local laws and restrictions.

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FCC / IC CLASS A DIGITAL DEVICE EMC VERIFICATION STATEMENT

NOTE: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC rules and Canadian ICES-003/NMB-003 regulation. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

Adam Equipment products have been tested with, and are always supplied with mains power adaptors which meet all legal requirements for the intended country or region of operation, including electrical safety, interference and energy efficiency. As we often update adaptor products to meet changing legislation it is not possible to refer to the exact model in this manual. Please contact us if you need specifications or safety information for your particular item. Do not attempt to connect or use an adaptor not supplied by us.

ADAM EQUIPMENT is an ISO 9001:2008 certified global company with more than 40 years' experience in the production and sale of electronic weighing equipment.

Adam products are predominantly designed for the Laboratory, Educational, Health and Fitness, Retail and Industrial Segments. The product range can be described as follows:

- Analytical and Precision Laboratory Balances
- Compact and Portable Balances
- High Capacity Balances
- Moisture analysers / balances
- Mechanical Scales
- Counting Scales
- Digital Weighing/Check-weighing Scales
- High performance Platform Scales
- Crane scales
- Mechanical and Digital Electronic Health and Fitness Scales
- Retail Scales for Price computing

For a complete listing of all Adam products visit our website at
www.adamequipment.com

Adam Equipment Co. Ltd. Maidstone Road, Kingston Milton Keynes MK10 0BD UK Phone: +44 (0)1908 274545 Fax: +44 (0)1908 641339 e-mail: sales@adamequipment.co.uk	Adam Equipment Inc. 1, Fox Hollow Rd. 06478 USA Phone: +1 203 790 4774 Fax: +1 203 792 3406 e-mail: sales@adamequipment.com	AE Adam GmbH. Instenkamp 4 D-24242 Felde Germany Phone +49 (0)4340 40300 0 Fax: +49 (0)4340 40300 20 e-mail: sales@adamequipment.co.za
Adam Equipment S.A. (Pty) Ltd. 7 Megawatt Road, Spartan EXT 22 Kempton Park, Johannesburg, Republic of South Africa Phone +27 (0)11 974 9745 Fax: +27 (0)11 392 2587 e-mail: sales@adamequipment.co.za	Adam Equipment (S.E. ASIA) PTY Ltd 2/71 Tacoma Circuit CANNING VALE 6155 Perth Western Australia Phone: +61 (0) 8 6461 6236 Fax +61 (0) 8 9456 4462 e-mail: sales@adamequipment.com.au	Adam Equipment (Wuhan) Co. Ltd. A Building East Jianhua Private Industrial Park Zhuanyang Avenue Wuhan Economic & Technological Development Zone 430056 Wuhan P.R.China Phone: + 86 (27) 59420391 Fax + 86 (27) 59420388 e-mail: info@adamequipment.com.cn

© Copyright by Adam Equipment Co. All rights reserved. No part of this publication may be reprinted or translated in any form or by any means without the prior permission of Adam Equipment.

Adam Equipment reserves the right to make changes to the technology, features, specifications and design of the equipment without notice.

All information contained within this publication is to the best of our knowledge timely, complete and accurate when issued. However, we are not responsible for misinterpretations which may result from the reading of this material.

The latest version of this publication can be found on our Website.

www.adamequipment.com

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© Copyright par Adam Equipment Co. Ltd. Tous droits réservés. Aucune partie de cette publication ne peut être reproduite ou traduite sous quelque forme ou par tout moyen, sans l'autorisation préalable d'Adam Equipment.
© copyright Adam Equipment Co. Ltd. Alle Rechte vorbehalten. Kein Teil dieser Publikation darf ohne vorherige Genehmigung von Adam Equipment in irgendeiner Form nachgedruckt oder übersetzt werden.
© Copyright by Adam Equipment Co. Ltd. Todos los derechos reservados. Ninguna parte de esta publicación puede ser reimpresa o traducida de ninguna forma o a través de ningún medio sin la previa autorización de Adam Equipment.
© Copyright di ADAM Equipment Ltd.