



Speed, Performance, Value - The Right Balance

- 9 weighing units
- Full range tare
- Backlit LCD display
- Stainless steel pan
- ShockProtect™ overload three-point protection
- Lock down mounting slot for Kensington® type lock and cable
- Dual tare keys
- Removable draught shield (120mm pan size only)
- Lightweight and portable
- Modern low profile design
- Colour coded, sealed keypad
- Low cost and simple to use
- External calibration
- Auto sleep / power down function to save battery life
- Auto power off
- AC Adapter
- Low battery indication
- Below balance weighing with hanger
- Weighing
- Below balance weighing



Technical Specifications

Maximum capacity	1750g
Readability	0.1g
Repeatability (S.D.)	0.1g
Linearity (\pm)	0.2g
Pan Size	120mm \emptyset

General Specifications

Model	CQT 1751GR
Units of measurement	Gram (g)
Power supply	12VAC 15mA
Draught shield dimensions	130 x 90mm
Net weight	0.9kg

Key Accessories

- In-Use Cover pack of 10
- Hard Carry Case with Lock
- Security Lock and Cable
- In-use Cover
- In-Use Cover pack of 5



Adam Equipment - offices worldwide

www.adamequipment.com