

Technical Specifications


Model	TBB 2610T	TBB 2610S
Maximum capacity	2610g	2610g
Readability	0.1g	0.1g
Tare range	+225g	-
Repeatability (S.D.)	0.2g	0.2g
Linearity (\pm)	0.2g	0.2g
Pan size	152mm \emptyset	

General Specifications

Units of measurement	Gram (g)
Calibration	External
Display	Graduated Beam
Housing	Aluminium
Overall dimensions (wxdxh)	505x110x160mm
Net weight	2.7kg

Key Accessories

- Security Lock and Cable
- Hard Carry Case with Lock



Adam Equipment - offices worldwide

www.adamequipment.com