

Technical Specifications

Maximum capacity	1750g
Readability	0.1g
Repeatability (S.D.)	0.1g
Linearity (\pm)	0.2g
Pan Size	120mm \emptyset

General Specifications

Model	CQT 1751GR
Units of measurement	Gram (g)
Power supply	12VAC 15mA
Draught shield dimensions	130 x 90mm
Net weight	0.9kg

Key Accessories

- In-Use Cover pack of 10
- Hard Carry Case with Lock
- Security Lock and Cable
- In-use Cover
- In-Use Cover pack of 5



Adam Equipment - offices worldwide

www.adamequipment.com