




MDW-160M Mechanical Person Scale



Speed, Performance, Value - The Right Balance

- Height rod as standard (measures up to 60-212 cm)
- Robust metal housing
- Anti-slip platform
- Adjustable zero point
- Wheels to allow easy movement
- 2 weighing units (kg, lb)
- Weighing
- Measuring height




www.adamequipment.com

Technical Specifications

Maximum capacity	160kg
Readability	100g
Repeatability (S.D.)	200g
Linearity (\pm)	200g
Platform size	375x275mm

General Specifications

Model	MDW-160M
Units of measurement	Kilogram (kg), Pound (lb)
Division of measurement	0.1cm
Housing	Metal
Overall dimensions (wxdxh)	530x275x1485mm
Net weight	15kg



Adam Equipment - offices worldwide

www.adamequipment.com