

MDW-160M Mechanical Person Scale



Technical Specifications

Maximum capacity	160kg
Readability	100g
Repeatability (S.D.)	200g
Linearity (\pm)	200g
Platform size	375x275mm

General Specifications

Model	MDW-160M
Units of measurement	Kilogram (kg), Pound (lb)
Division of measurement	0.1cm
Housing	Metal
Overall dimensions (wxdxh)	530x275x1485mm
Net weight	15kg



Adam Equipment - offices worldwide

www.adamequipment.com