


MDW-250L Person Scale



Speed, Performance, Value - The Right Balance

- Height rod as standard (measures up to 71-210 cm)
- Robust metal housing
- Anti-slip platform
- Large red LED display
- Colour coded, sealed keypad
- Hold function
- External calibration
- Full range tare
- Overload protection
- AC Adapter
- Rechargeable battery / mains adapter supplied as standard
- 2 weighing units (kg, lb)
- Precision load cell technology
- Simple user-friendly operation
- Weighing
- Measuring height



Technical Specifications

Maximum capacity	250kg
Readability	100g
Repeatability (S.D.)	0.2kg
Linearity (\pm)	200g
Platform size	375x275mm

General Specifications

Model	MDW-250L
Tare functions	Yes
Units of measurement	Kilogram (kg), Pound (lb)
Functions	Hold and Tare functions
Calibration	External
Display	LCD
Power supply	Internal re-chargeable 6V battery / 220 VAC through adapter
Operating temperature	5°C - 40°C
Height range	71cm - 210cm
Division of measurement	0.5cm
Other features	Automatic backlight, overload protection
Housing	Metal, Indicator: ABS Plastic
Overall dimensions (wxdxh)	535x275x940mm
Net weight	13kg



Adam Equipment - offices worldwide

www.adamequipment.com