

Technical Specifications

Maximum capacity	250kg
Readability	100g
Repeatability (S.D.)	0.2kg
Linearity (\pm)	200g
Platform size	375x275mm

General Specifications

Model	MDW-250L
Tare functions	Yes
Units of measurement	Kilogram (kg), Pound (lb)
Functions	Hold and Tare functions
Calibration	External
Display	LCD
Power supply	Internal re-chargeable 6V battery / 220 VAC through adapter
Operating temperature	5°C - 40°C
Height range	71cm - 210cm
Division of measurement	0.5cm
Other features	Automatic backlight, overload protection
Housing	Metal, Indicator: ABS Plastic
Overall dimensions (wxdxh)	535x275x940mm
Net weight	13kg



Adam Equipment - offices worldwide

www.adamequipment.com