

Technical Specifications


| Model | CBK 4 | CBK 8H | CBK 8 | CBK 16 | CBK 32 | CBK 48 |
|----------------------|-----------|--------|-------|--------|--------|--------|
| Maximum capacity | 4000g | 8000g | 8000g | 16kg | 32kg | 48kg |
| Readability | 0.1g | 0.1g | 0.2g | 0.5g | 1g | 2g |
| Repeatability (S.D.) | 0.1g | 0.1g | 0.2g | 0.5g | 1g | 2g |
| Linearity (±) | 0.2g | 0.2g | 0.4g | 1g | 2g | 4g |
| Platform | 225x275mm | | | | | |

General Specifications

| | |
|----------------------------|--|
| Units of measurement | Gram (g), Kilogram (kg) |
| Interface | RS-232 bi-directional Interface |
| Calibration | Automatic external |
| Display | 6 digits LCD displays |
| Battery | Internal Rechargeable battery up to 90 hours |
| Power supply | 230VAC 50/60Hz 120VAC available |
| Operating temperature | -10°C 40°C |
| Real-time clock | Real-time clock |
| Housing | ABS Plastic |
| Overall dimensions (wxdxh) | 315x355x110mm |
| Net weight | 4.1kg |

Key Accessories

- Adam Printer
- Thermal Paper Roll for Adam Printer
- In-use Cover
- In-Use Cover pack of 10
- Weigh Below Hanger (factory fitted)
- RS-232 to PC Cable
- RS-232 to USB Interface Cable (must be ordered with 3074010266 RS232 cable)
- AdamDU - Data Collection Program
- Hard Carry Case with Lock
- In-Use Cover pack of 5
- ATP Thermal Printer
- Thermal Paper for ATP



Adam Equipment - offices worldwide

www.adamequipment.com